



Botico

Making products, delivering Success

BOTSYN MEL

CHARACTERISTICS

Appearance:

Pale white powder

Base:

Melamine

Solids:

93 (+/- 2) %

pH (@ 10% aqueous soln.):

5.5 - 7.5

Charge:

Anionic

Light Fastness:

Medium

REACH COMPLIANT

Date of issue

Melamine based resin for retanning.

BOTSYN MEL is suitable for retanning leathers where it gives full mellow leathers.

BOTSYN MEL helps in filling loose structured areas of leather

It improves buffing property of Nubuck and suede leathers.

Application

The recommended quantities of BOTSYN MEL are 3 - 5 % depending on article being processed.

02-08-2023

Shelf life

The product can be stored in its original packaging between 5°C and 35°C for 12 months.

The information and recommendations given are provided on the basis of our best current knowledge, supported by applied experience and accurate testing. However, this information does not constitute a guarantee or responsibility for product use, which must always be checked and evaluated by each individual client.